



2000W USB 2.0 Port Skin Analysis Machine 7200K Skin Analyzer Device With Magnifier

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: Bitmoji
- Certification: ISO CE
- Model Number: A6
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



Product Specification

- Name: Skin Analysis Machine
- Packing Size: L68xW49xH51cm
- Color Temperature: 7200K
- Magnifier: Yes
- Camera Sensor Model: CCD1/1.7
- Weight: 2.5kg
- Detect Problem: 12 Skin Problem Detected
- Type: Portable
- Megapixel: 2000W
- Highlight: 2000W Skin Analysis Machine, 2.0 Port Skin Analysis Machine, 7200K Skin Analyzer Device



More Images





Product Description

2000W USB 2.0 Port Skin Analyzer Facial Skin Analyzer With 7200K Color Temperature And Magnifier

Product Description:

The machine is compact and has a packing size of L68xW49xH51cm, which makes it easy to carry and store. The handheld skin analyzer is perfect for beauty salons, dermatology clinics, and personal use.

The Face Analyzer Machine operates on an RGB+UV operating system, which provides a clear and accurate image of the skin. The machine is equipped with a magnifier that helps you observe the skin texture, pores, pigmentation, and other details that are not visible to the naked eye. The magnifier is easy to adjust, and you can customize the settings to suit your needs.

The Facial Dermabrasion Machine is compatible with Windows XP/Vista/7/8/10 operating systems, which makes it easy to integrate with your computer. You can view the skin analysis report on your computer screen and save it for future reference. The machine is easy to use and does not require any special training or skills.

The Handheld Skin Analyzer is a versatile and user-friendly device that helps you understand your skin type and condition. The machine is suitable for all skin types and can be used on the face, neck, and other body parts. The machine is ideal for people who want to take control of their skin health and beauty.

In conclusion, the Facial Dermabrasion Machine is an excellent investment for anyone who wants to analyze their skin condition accurately and effectively. The machine is easy to use, reliable, and provides accurate results. With its advanced features and user-friendly design, the Handheld Skin Analyzer is a must-have device for all beauty enthusiasts and professionals.

Features:

Product Name: Skin Analysis Machine

Operating System: RGB +UV

Computer Requirement: Windows XP/Vista/7/8/10

Working Environment: 10 – 35 Celsius Degrees

Usb: 2.0 Port

Magnifier Rate: 2 X

Other names for this product include: Body Analyzer Machine, Digital Skin Analyzer, Digital Skin Camera Microscope.

Technical Parameters:

Type	Portable
Computer Requirement	Windows XP/Vista/7/8/10
Pdf Print	YES
Power Supply	AC 100-240V, 50/60Hz
Color Temperature	7200 K
Operating System	RGB +UV
Magnifier Rate	2 X
Magnifier	Yes
Detect Problem	12 Skin Problem Detected
Packing Size	L68xW49xH51cm

A6-AI Intelligent Image Instrument Product Specifications

Your skin diagnostic expert! Professional Skin Care Management Solution



Device Name	A6-AI intelligent skin tester
Model	ZMLH-A6
Pixel	Industrial grade 36 Million pixel
Spectrum	RGB white light; Positive polarized light; Negative polarized light; Wood's light; UV light; Red light; Brown light; Mixed light
Shading Method	Semi-open
Product Material	Industrial Grade ABS
Power Supply	AC 100-260V, 50/60Hz
Product Dimensions	Unfolded: 636 x 521 x 489mm; Folded: 361 x 521 x 489mm
Central Processor	Rockchip RK3568 ; Quad-core Cortex-A55
Motherboard	R10-S6810 Pmotherboard, running Android 11
Main Memory	Dual-channel LP DD R4, 4G
Storage	MMC 5.1, 32G
Operating Method	multi-point capacitive touch
Screen Features	Foldable 90 degrees
WIFI	Built-in dual-band WIFI (2.4G, 5G)
HDMI	1 port (optional)
USB	2 ports
Screen Ratio	16:9
Screen Size	15.6 inches
Screen Resolution	1920*1080
Net Weight	9.5 KG
Gross Weight	14.6 KG
Packaging Dimensions	640 x 570 x 605mm
Packaging Material	Corrugated paper + pearl cotton
Other Accessories	Handheld skin analyzer, power cord

Bitmoji®

A6-AI intelligent skin tester



It currently supports **16** languages

**More language
versions are
coming soon!**



Look forward to **more languages!**





Applications:

One of the main benefits of this product is its ability to accurately assess the condition of the skin. It can detect issues such as dryness, oiliness, wrinkles, and blemishes. This allows for the creation of personalized skincare routines tailored to individual needs.

The Facial Skin Analyzer Machine is also useful in providing visual evidence of the effectiveness of skincare treatments. By using the magnifier rate of 2X, users can see the results of procedures such as facial dermabrasion machine and chemical peels.

Thanks to its place of origin in China Guangzhou, this product benefits from the expertise of Chinese manufacturers in the field of skincare technology. It is made with high-quality materials and is designed to withstand heavy use.

The product is also compatible with a wide range of computer systems, including Windows XP/Vista/7/8/10. This allows for easy integration into existing technology setups.

Finally, the color temperature of 7200 K ensures that the images captured by the Facial Skin Analyzer Machine are accurate and true to life. Users can rely on the product to provide them with an in-depth analysis of their skin condition, and make informed decisions about their skincare routine.

Overall, the Facial Skin Analyzer Machine is a versatile and powerful tool for anyone looking to take their skincare to the next level. Its applications are wide-ranging, from personal use at home to professional use in dermatology clinics. With the ability to accurately assess and track skin condition, it is an invaluable asset for anyone looking to achieve healthy, glowing skin.

Pdf Print: YES

Magnifier Rate: 2 X

Power Supply: AC 100-240V, 50/60Hz

Computer Requirement: Windows XP/Vista/7/8/10

Color Temperature: 7200 K

Customization:

Our Skin Analysis Machine is a top-of-the-line dermatological equipment that provides accurate skin analysis and diagnosis to help you determine the best facial skin care regimen. It is made in China Guangzhou, ensuring the highest quality of materials and craftsmanship. With a weight of only 2.5kg, it is lightweight and portable, making it convenient for use in various locations. The machine has a color temperature of 7200 K, providing natural and accurate lighting for better skin analysis results.

Equipped with an RGB + UV operating system, this Acne Removal Machine can detect various skin issues such as pigmentation, acne, and wrinkles. It also has a PDF print function, allowing you to easily print out the analysis report for your clients or personal records.

The working environment of the machine should be between 10 – 35 Celsius Degrees to ensure optimal performance. Our Skin Analysis Machine is a must-have for any professional in the beauty industry who wants to provide the best possible facial skin care services.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou

